

## Enteric Pathogen Panel



<b>Alternate Name:</b>	Epic: LAB311025, MLAB: EPV
<b>Performing Lab:</b>	New Hanover
<b>Test Includes:</b>	Campylobacter group, Salmonella species, Shigella species, Vibrio group, Yersinia enterocolitica, Norovirus GI/GII and Rotavirus A)
<b>Specimen Container:</b>	Screw cap container
<b>Minimum Volume Required:</b>	2.0 grams
<b>Testing Availability</b>	<b>Routine:</b> 24 hours/day <b>Stat:</b> No
<b>Turnaround Time:</b>	6 – 12 hours
<b>Special Handling:</b>	Specimens must be brought to the laboratory within: 4 hours- Ambient specimens 24 hours- Refrigerated specimens 48 hours- specimens in Cary Blair
<b>Patient Preparation:</b>	None
<b>Specimen Stability:</b>	N/A
<b>Reference Range:</b>	N/A
<b>Critical Value:</b>	N/A
<b>CPT Code:</b>	87045
<b>Testing Methodology:</b>	Amplified Molecular
<b>Causes for Specimen Rejection:</b>	Improper labeling Specimen container must not be contaminated on the outside. Specimens must be delivered to the laboratory within the above specimen handling time limits.
<b>Other Comments:</b>	N/A
<b>Clinical Significance:</b>	N/A

