

Acceptable Tube Types

Hepatitis C Virus Antibody (IgG)



Alternate Name:	Anti-HCV, Hepatitis Non-A, Non-B
Performing Lab:	New Hanover
Specimen Container:	Pink or Lavender top tube(EDTA), Yellow serum separator tube, or a red top tube.
Minimum Volume Required:	1.0 mL
Testing Availability	Routine: 24 hours/day Stat: No
Turnaround Time:	Routine: 4 hours
Special Handling:	None
Patient Preparation:	None
Specimen Stability:	7 Days at 2-8° C
Reference Range:	Nonreactive
Critical Value:	None
CPT Code:	86803
EPIC Test Number:	LAB868
Testing Methodology:	Chemiluminescence
Causes for Specimen Rejection:	Improper labeling
Other Comments:	
Clinical Significance:	Assess exposure to hepatitis C virus infection.