

Test Name:	Stix Test, Urine
Performing Lab:	New Hanover and Cape Fear
Clinical Significance:	
Tests Included:	pH, Blood, Ketone, and Protein (When ordering STIX testing, a maximum of three of the four tests may be selected)
Reference Range:	pH: 4.5 – 8.0 Blood: Negative Ketone: Negative Protein: Trace
Critical Value:	
Specimen Container:	Clean screw-cap container
Minimum Volume:	1 mL
Special Handling:	None
Preparation of Patient:	None
Specimen Stability:	
Causes for Specimen Rejection:	Improper labeling
Testing Methodology:	Reagent strips
Test Performed:	Daily, 24 hours
Turnaround Time:	Routine: 4 hours Stat: < 1 hour
CPT Code:	