

Stool for WBC



Alternate Name:	Fecal Lactoferrin
Performing Lab:	New Hanover
Specimen Container:	Clean screw capped container
Minimum Volume Required:	1 mL
Testing Availability	Routine: 24 hours / day Stat: No
Turnaround Time:	24 hours
Special Handling:	None
Patient Preparation:	None
Specimen Stability:	N/A
Reference Range:	Negative for WBC
Critical Value:	N/A
CPT Code:	N/A
Testing Methodology:	Immunochromatographic
Causes for Specimen Rejection:	Improper labeling. Purgatives other than saline or Fleet.
Other Comments:	None
Clinical Significance:	Conditions associated with fecal leukocytes include diffuse antibiotic associated colitis, ulcerative colitis, shigellosis, salmonellosis, <i>Campylobacter</i> , and <i>Yersinia</i> infection.