



Calcium, Urine

Alternate Name:	None	
Performing Lab:	New Hanover	
Specimen Container:	Random:	Urine Specimen Cup
	24 Hour:	24 Hour Urine Collection Container
Minimum Volume Required:	Random:	1.0 mL
	24 Hour:	Not Applicable
Testing Availability	Routine:	24 hours / day
	Stat:	No
Turnaround Time:	4 hours	
Special Handling:	24-Hour collection: Requisition must state date and time collection started and date and time collection finished. Urine should be refrigerated during collection. Do not use preservatives.	
Patient Preparation:	None	
Specimen Stability:		
Reference Range:	Random:	None
	24 Hour:	42 – 353 mg/24 hr
Critical Value:	N/A	
CPT Code:	Random Urine:	82310
	24 Hour Urine:	82340
Testing Methodology:	Two-point rate	
Causes for Specimen Rejection:	Improper labeling Use of preservatives	
Clinical Significance:	N/A	